

SUPPORTING SCIENTISTS DEVELOPING MOLECULAR DIAGNOSTICS

NGS- AND AMPLIFICATION-BASED SOLUTIONS



IT'S TIME TO **THINK DIFFERENTLY**

The COVID-19 pandemic elevated many of the challenges faced by clinical labs, requiring new and innovative solutions to address them. Technology development is happening faster than ever before, increasing the need for innovation and thinking differently about how diagnostics should be developed, manufactured and deployed.

Many scientists know NEB as a trusted reagent provider to the life science community. What many do not know is that we also offer products and support for scientists who are working to develop diagnostics. Extensive molecular biology and enzymology experience enables us to help customers solve the challenges inherent in technology development and ultimately in scale-up and commercialization. For scientists looking for long term, ambient-stored products, NEB's diverse portfolio of lyophilized and lyo-compatible products, as well as custom lyophilization capabilities are also available.

LEARN MORE AT

www.neb.com/MDx

DEVELOPING A **LIFE SCIENCE PRODUCT OR TECHNOLOGY?**

Contact our Customized Solutions Team to build your personalized solution.

With 50 years of experience in enzymology and reagent manufacturing, NEB's Customized Solutions Team will help you access novel products, meet quality specifications, speed time to market, and streamline your supply chain, allowing you to focus more on what matters most — innovation.

AREAS WE CAN HELP YOU:



Bulk Packaging



Custom Manufacturing



Packaging & Private Label



GMP-grade* Reagents



Lyophilization Capabilities

GET STARTED AT

www.neb.com/customizedsolutions

GMP-grade is a branding term NEB uses to describe reagents manufactured at NEB's Rowley facility. The Rowley facility was designed to manufacture reagents under more rigorous infrastructure and process controls to achieve more stringent product specifications and customer requirements. Reagents manufactured at NEB's Rowley facility are manufactured in compliance with ISO 9001 and ISO 13485 quality management system standards. However, at this time, NEB does not manufacture or sell products known as Active Pharmaceutical Ingredients (APIs), nor does NEB manufacture its products in compliance with all of the Current Good Manufacturing Practice regulations.



NGS-BASED APPLICATIONS

Sequencing is enabling scientists to make rapid advances in epidemiology and surveillance, basic and clinical research, and diagnostics. A fast-growing number of methods are being developed to address whole genome, as well as targeted approaches. In all cases, streamlined workflows that result in high-quality, high yield libraries are critical towards optimizing your next generation sequencing (NGS) results. The table below summarizes some of the NGS products available for scientists who are working to develop diagnostics.

APPLICATION	PRODUCTS		NOTES
DNA Sequencing	DNA library preparation with enzymatic DNA fragmentation	NEBNext UltraExpress® FS DNA Library Prep Kit (NEB #E3340) NEBNext Ultra™ II FS DNA Library Prep Kit for Illumina® (NEB #E7805) NEBNext Ultra II FS DNA PCR-free Library Prep Kit for Illumina (NEB #E7730) NEBNext Ultra II FS DNA Module (NEB #E7810) NEBNext UltraShear® FFPE DNA Library Prep (NEB #E6655) NEBNext Ultra II Q5® Master Mix (NEB #M0544) NEBNext® Adaptors and Primers (www.neb.com/oligos)	<ul style="list-style-type: none"> • Range of packaging options including 96-well plates • Component and kit customization • Kitting • Automation-compatible • Streamlined workflows
	DNA library preparation for cfDNA or pre-sheared DNA	NEBNext UltraExpress DNA Library Prep Kit (NEB #E3325) NEBNext Ultra II DNA Library Prep Kit for Illumina (NEB #E7645) NEBNext Ultra II DNA PCR-free Library Prep Kit for Illumina (NEB #E7410) NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546) NEBNext Ultra II Ligation Module (NEB # E7595) NEBNext Ultra II Q5 Master Mix (NEB #M0544) NEBNext Adaptors and Primers (www.neb.com/oligos)	
	DNA methylation sequencing	NEBNext Enzymatic Methyl-seq Kit (NEB #E7120) NEBNext Enzymatic 5hmC-seq Kit (NEB #E3350)	
NGS-BASED APPLICATIONS	Depletion of abundant RNAs and enrichment of poly(A)-tailed mRNA	NEBNext Poly(A) mRNA Magnetic Isolation Module (NEB #E7490) NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat) (NEB #E7400) NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) (NEB #E7750) NEBNext rRNA Depletion Kit (Bacteria) (NEB #E7850) NEBNext RNA Depletion Core Reagent Set (for customized depletion) (NEB #E7865)	<ul style="list-style-type: none"> • Range of packaging options including 96-well plates • Component and kit customization • Kitting • Automation-compatible • Streamlined workflows
	RNA library preparation	NEBNext UltraExpress™ RNA Library Prep Kit (NEB #E3330) NEBNext Ultra II Directional RNA Library Prep Kit for Illumina (NEB #E7760) NEBNext Ultra II RNA First Strand Synthesis Module (NEB #E7771) NEBNext Ultra II Directional RNA Second Strand Synthesis Module (NEB #E7550) NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546) NEBNext Ultra II Ligation Module (NEB #E7595) NEBNext Ultra II Q5 Master Mix (NEB #M0544) NEBNext Adaptors and Primers (www.neb.com/oligos) NEBNext Immune Sequencing Kit (Human) (NEB #E6320)	
	Virus sequencing and detection	NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina) (NEB #E7658) NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina) (NEB #E7650) NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies) (NEB #E7660) NEBNext ARTIC SARS-CoV-2 RT-PCR Module (NEB #E7626) Q5 Hot Start High-Fidelity 2X Master Mix (NEB #M0494) NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546) Blunt/TA Ligase Master Mix (NEB #M0367) NEBNext Quick Ligation Module (NEB #E6056) NEBNext Ultra II RNA Library Prep Kit for Illumina (NEB #E7770) NEBNext Adaptors and Primers (www.neb.com/oligos)	



AMPLIFICATION-BASED APPLICATIONS

As the first company to sell *Taq* DNA Polymerase, NEB has a long history in the development of reliable and convenient tools for amplification, and offers a large selection of products for PCR, qPCR, RT-qPCR and isothermal amplification. Our expertise in this area has allowed us to develop optimized enzymes for a variety of applications, that are available for scientists working to develop diagnostics. The table below summarizes some of the products available.

APPLICATION	PRODUCTS	PRODUCT NOTES	CUSTOM FORMULATIONS AVAILABLE
PCR APPLICATIONS	DNA, Dye Luna® Universal qPCR Master Mix (NEB #M3003) DNA, Probe Luna Universal Probe qPCR Master Mix (NEB #M3004)	<ul style="list-style-type: none"> Compatible with automation and reaction miniaturization Room temperature stable for ≥ 24 hours 	<ul style="list-style-type: none"> Blue-dye-free Lyo-compatible
	RNA (1-step), Dye Luna Universal One-Step RT-qPCR Kit (NEB #E3005) RNA (1-step), Probe Luna Probe One-Step RT-qPCR 4X Mix with UDG (NEB #M3019) Luna Probe One-Step RT-qPCR 4X Mix with UDG (No ROX) (NEB #M3029) Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006) Luna Probe One-Step RT-qPCR Kit (No ROX) (NEB #E3007) Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit (NEB #E3019) LyoPrime Luna™ Probe One-Step RT-qPCR Mix with UDG (NEB #L4001)	<ul style="list-style-type: none"> Luna WarmStart RT paired with Hot Start <i>Taq</i> increases reaction specificity and robustness Compatible with automation and reaction miniaturization Room temperature stable for ≥ 24 hours Lyophilized format (NEB #L4001) 	<ul style="list-style-type: none"> Blue-dye-free Lyo-compatible Add primers and probes (NEB #L4001)
	RNA (2-step) LunaScript® RT SuperMix Kit (NEB #E3010/M3010)	<ul style="list-style-type: none"> Novel engineered RT Single-tube format 13-minute cDNA synthesis protocol 	<ul style="list-style-type: none"> Blue-dye-free Primer-free
	Master Mixes Q5® Hot Start High-Fidelity 2X Master Mix (NEB #M0494) Q5 High-Fidelity 2X Master Mix (NEB #M0492) Standalone Enzyme & Buffer Q5 Hot Start High-Fidelity DNA Polymerase (NEB #M0493) Q5 High-Fidelity DNA Polymerase (NEB #M0491) Q5 Blood Direct 2X Master Mix (NEB #M0500)	<ul style="list-style-type: none"> ~280X fidelity of <i>Taq</i> Consistent, fast, reliable performance Compatible with automation and reaction miniaturization Room temperature stable for ≥ 24 hours 	<ul style="list-style-type: none"> High conc. Glycerol-free Custom mixes
	Hot Start <i>Taq</i> DNA Polymerase (NEB #M0495) Hot Start <i>Taq</i> 2X Master Mix (NEB #M0496)	<ul style="list-style-type: none"> Unique aptamer-based enzyme control supports fast protocols 	<ul style="list-style-type: none"> High conc. Glycerol-free
	WarmStart® Colorimetric LAMP 2X Master Mix (DNA & RNA) (NEB #M1800) WarmStart Colorimetric LAMP 2X Master Mix with UDG (NEB #M1804)	<ul style="list-style-type: none"> Fast, clear pink-to-yellow visible detection of amplification Results in approximately 30 minutes Automation-compatible when coupled with absorbance plate reader 	
	SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit (NEB #E2019)	<ul style="list-style-type: none"> Simple, colorimetric detection of amplification of SARS-CoV-2 nucleic acid 	
	WarmStart LAMP Kit (DNA & RNA) (NEB #E1700) WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG) (NEB #M1708) WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG) (NEB #E1708) LyoPrime WarmStart Fluorescent LAMP/RT-LAMP Mix with UDG (NEB #L4401)	<ul style="list-style-type: none"> Master mix for LAMP and RT-LAMP workflows Supports multiple detection methods, including fluorescence and turbidity Lyophilized format (NEB #L4401) 	<ul style="list-style-type: none"> Lyo-compatible High conc. Add primers (NEB #L4401)
	<i>Tte</i> UvrD Helicase (NEB #M1202)	<ul style="list-style-type: none"> Improves specifically of fluorescent LAMP reactions 	<ul style="list-style-type: none"> High conc.
	<i>Bst</i> 2.0 WarmStart DNA Polymerase (NEB #M0538) <i>Bst</i> 2.0 WarmStart DNA Polymerase (Glycerol-free) (NEB #M0402) <i>Bst</i> 2.0 DNA Polymerase (NEB #M0537)	<ul style="list-style-type: none"> Improved reaction properties compared to wild-type <i>Bst</i> DNA Polymerase Increased dUTP tolerance enables carryover prevention 	<ul style="list-style-type: none"> Glycerol-free (NEB #M0402) High conc.



APPLICATION	PRODUCTS	PRODUCT NOTES	CUSTOM FORMULATIONS AVAILABLE	
ISOTHERMAL APPLICATIONS	LAMP/ RT-LAMP	Bst 3.0 DNA Polymerase (NEB #M0374)	<ul style="list-style-type: none"> DNA binding domain fusion supports robust performance Significantly increased RT activity up to 72°C enables single enzyme RT-LAMP 	<ul style="list-style-type: none"> Glycerol-free High conc.
		WarmStart RTx Reverse Transcriptase (NEB #M0380) WarmStart RTx Reverse Transcriptase (Glycerol-free) (NEB #M0439)	<ul style="list-style-type: none"> <i>In-silico</i> designed RT for RT-LAMP with reversibly-bound aptamer that inhibits activity below 40°C 	<ul style="list-style-type: none"> Glycerol-free (NEB #M0439) High conc.
	Rolling Circle Amplification/ Whole Genome Amplification	phi29-XT RCA Kit (NEB #E1603) phi29-XT WGA Kit (NEB #E1604)	<ul style="list-style-type: none"> Includes phi29-XT DNA polymerase, an engineered polymerase that generates high product yield in a short time with improved thermostability and sensitivity 	
	Strand Displacement	Nt.BstNBI (NEB #R0607) WarmStart Nt. BstNBI (NEB #R0725) WarmStart Afu Uracil-DNA Glycosylase (UDG) (NEB #M1282)	<ul style="list-style-type: none"> High purity, high quality nicking endonuclease Warmstart Afu UDG digests uracil-containing amplification at temps over 40°C 	<ul style="list-style-type: none"> Glycerol-free High conc.
	Helicase-dependent Amplification	IsoAmp® II Universal tHDA Kit (NEB #H0110)	<ul style="list-style-type: none"> Requires only two primers Produces short, discrete DNA products 	
	Recombinase Polymerase Amplification	T4 UvsX DNA Recombinase (Glycerol-free) (NEB #M3081) T4 UvsY Protein (Glycerol-free) (NEB #M3082) Bsu DNA Polymerase Large Fragment (NEB #M0330) T4 Gene 32 Protein (NEB #M0300)	<ul style="list-style-type: none"> Core components of RPA reaction Can increase yield and efficiency of amplification reactions 	<ul style="list-style-type: none"> High conc. Glycerol-free
	Other	Deoxynucleotide (dNTP) Solution Mix (NEB #N0447)		<ul style="list-style-type: none"> Custom conc.
		Nuclease-free water (NEB #B1500)		
		Antarctic Thermolabile UDG (NEB #M0372)	<ul style="list-style-type: none"> Unique thermolabile version is completely inactivated in typical isothermal and RT-qPCR workflows 	<ul style="list-style-type: none"> High conc.
		Proteinase K, Molecular Biology Grade (NEB #P8107)		
	Thermolabile Proteinase K (NEB #P8111)	<ul style="list-style-type: none"> Unique thermolabile version is completely inactivated in typical isothermal and RT-qPCR workflows 	<ul style="list-style-type: none"> Custom conc. 	

To inquire about bulk and/or custom formats, or GMP-grade* reagent production, contact us at

www.neb.com/customizedsolutions

USA

New England Biolabs, Inc.
Telephone (978) 927-5054
Toll Free (USA Orders) 1-800-632-5227
Toll Free (USA Tech) 1-800-632-7799
info@neb.com

Australia & New Zealand

New England Biolabs (Australia) PTY
Telephone: 1800 934 218 (AU)
info.au@neb.com
Telephone: 0800 437 209 (NZ)
info.nz@neb.com

Canada

New England Biolabs, Ltd.
Toll Free: 1-800-387-1095
info.ca@neb.com

China

New England Biolabs (Beijing), Ltd.
Telephone: 010-82378265/82378266
info@neb-china.com

France

New England Biolabs France
Telephone: 0800 100 632
info.fr@neb.com

Germany & Austria

New England Biolabs GmbH
Free Call: 0800/246 5227 (Germany)
Free Call: 00800/246 52277 (Austria)
info.de@neb.com

Japan

New England Biolabs Japan, Inc.
Telephone: +81 (0)3 5669 6191
info.jp@neb.com

Singapore

New England Biolabs, Pte. Ltd.
Telephone: +65 638 59623
sales.sg@neb.com

Republic of Korea

New England Biolabs Korea Ltd.
Toll Free: +82 (70) 47318478
info.kr@neb.com

United Kingdom

New England Biolabs (UK), Ltd.
Call Free: 0800 318486
info.uk@neb.com



Products and content are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. (NEB). The use of trademark symbols does not necessarily indicate that the name is trademarked in the country where it is being read; it indicates where the content was originally developed. See www.neb.com/trademarks. The use of these products may require you to obtain additional third-party intellectual property rights for certain applications. For more information, please email busdev@neb.com.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

© Copyright 2024, New England Biolabs, Inc.; all rights reserved.

